



INTERX 1000

Personal Device

Product# 1243-001 / SIN A 72



FEATURES

- **Displays** frequency setting and intensity settings
- **5 settings** include chronic and acute pre-sets
- **Lightweight, portable** with 2 AA DC Alkaline Batteries power source
- **InterX non-invasive interactive neurostimulation (NIN) technology** listens to biofeedback in order to responsively deliver an interactive and painless waveform that is pulsed, damped, and bi-phasic sinusoidal

SPECIFICATIONS

- **Size:** 215mm x 52mm x 43mm
- **Weight:** 6.5 ounces (185 grams)
- **Frequency:** 15-480 Hz

MKTG-1000-PRDCRD

InterX® Electrode Accessories

InterX® Electrode Accessories are a critical part of the InterX® non-invasive interactive neurostimulation (NIN) system. Electrode Accessories can be connected with InterX® Devices to deliver enhanced treatment as well as dynamic and unattended treatment. Each of these Electrode Accessories is compatible with the InterX® 1000 Personal Device.



3x3 DUAL FLEX ARRAY

Product# 1376-001 / SIN A 72

Use for neck, spine, joints, and extremities, and for unattended treatment.



4x4 FLEX ARRAY

Product# 1299-001 / SIN A 72

Use for spine, hamstring, and quad, and for unattended treatment.



CLASSIC

Product# 1569-001 / SIN A 72

Use for scanning face, hand, wrist, toes, ankles. Made for small treatment areas, giving better skin contact.



COMB

Product# 1070-001 / SIN A 72

Use for scalp and where hair may interfere with electrode contact.



COSMETIC

Product# 1075-001 / SIN A 72

Use for facial treatment.



DOME

Product# 1295-001 / SIN A 72

For large surface areas, concentrated treatment for dynamic application and unattended use.



SOFT TISSUE — Universal

Product# 1674-001 / SIN A 72

Use for soft tissue manipulation and massage for fascia restrictions and scar management.

See treatment guides @ www.interx.com