

INTERX®

INTERX 5002

Pro Device

Product# 1146-001 / SIN A 72



FEATURES

- **Displays skin impedance readings, skin contact verification and timer** during active treatment allowing user to identify optimal treatment points and treatment progress
- **Display** screen navigation of settings, intensity, battery life, and volume and brightness options
- **14 settings**—chronic and acute frequencies and pre-sets
- **Lightweight, portable** with 4 AA DC Alkaline Batteries power source
- **InterX non-invasive interactive neurostimulation (NIN) technology** listens to biofeedback in order to responsively deliver an interactive and painless waveform that is pulsed, damped, and bi-phasic sinusoidal

SPECIFICATIONS

- **Size:** 220mm x 63mm x 45mm
- **Weight:** Approx. 11.6 oz (330 g)
- **Frequency:** 15-480 Hz

MKTG-5002-PRDCRD

InterX® Electrode Accessories

InterX® Electrode Accessories are a critical part of the InterX® non-invasive interactive neurostimulation (NIN) system. Electrode accessories can be connected with InterX® Devices to Electrode accessories to deliver enhanced treatment as well as dynamic and unattended treatment. All Electrode Accessories are compatible with the InterX 5002.



3x3 DUAL FLEX ARRAY

Product# 1376-001 / SIN A 72

Use for neck, spine, joints, and extremities, and for unattended treatment.



4x4 FLEX ARRAY

Product# 1299-001 / SIN A 72

Use for spine, hamstring, and quad, and for unattended treatment.



CLASSIC

Product# 1569-001 / SIN A 72

Use for scanning face, hand, wrist, toes, ankles. Made for small treatment areas, giving better skin contact.



COMB

Product# 1070-001 / SIN A 72

Use for scalp and where hair may interfere with electrode contact.



COSMETIC

Product# 1075-001 / SIN A 72

Use for facial treatment.



DOME

Product# 1295-001 / SIN A 72

For large surface areas, concentrated treatment for dynamic application and unattended use.



SOFT TISSUE — Universal

Product# 1674-001 / SIN A 72

Use for soft tissue manipulation and massage for fascia restrictions and scar management.



SOFT TISSUE — Small

Product# 1538-001 / SIN A 72

Use for hands and feet. For use with InterX 5002 Pro device only.

See treatment guides @ www.interx.com